

LOABeads Protein A – Quick Protocol

For all separation steps, use the LOABeads MagSep Cube or LOABeads MagSep 15/50 magnets. Only remove liquid when beads are fully separated by magnet.

Purification of 1–2 mg IgG from serum or ascites

1. Dispense 500 μ l well-mixed bead suspension into a test tube.
2. Wash 1x with 500 μ l PBS.
3. Add 500 μ l PBS + 250–500 μ l serum or ascites.
4. Mix for 30–60 min.
5. Wash at least 3x with 500 μ l PBS.
6. Add 300–500 μ l 60 mM citrate, pH 3.0.
7. Mix for 1 min to elute.
8. Add 30–50 μ l 2 M Tris-HCl, pH 9.0, to the saved supernatant for neutralization.
9. Elute a 2nd time if needed.
10. Regenerate beads and store in PBS, 20% EtOH.

Depletion of IgG from serum

1. Perform steps 1–4 above.
2. Save the supernatant, which represents the IgG-depleted serum.
3. Regenerate beads and store in PBS, 20% EtOH.

Immunoprecipitations

- Use 50–200 μ l 10% LOABeads Protein A suspension (5–20 μ l settled beads) per immunoprecipitation.
- Perform all steps throughout the immunoprecipitation as with classical, non-magnetic protein A agarose, except use the LOABeads MagSep Cube magnet for separations.