



Model: CCL-170/240_-_-HHS

CelCulture®

CO₂ Incubators with High Heat SterilizationCultivating a Culture of Safety and Efficiency





Welcome to Esco

Esco's Vision is to provide enabling technologies for scientific discoveries to make human lives healthier and safer.

The Esco Lifesciences Group is committed to deliver innovative solutions for the clinical, life sciences, research, industrial, laboratory, pharmaceutical, and IVF community. With the most extensive product line in the industry, Esco have passed a number of international standards and certifications. Esco represents innovation and forward-thinking designs, that are of the highest standard quality since 1978.

Availability and Accessibility. Esco has headquarters in Singapore, Indonesia, and Philippines, with manufacturing facilities are located in Asia and Europe. Research and Development (R&D) is conducted worldwide spanning the US, Europe and Asia. Sales, services and marketing subsidiaries are located in 42 major markets including US, UK, Japan, China and India. Esco regional distribution centers are located in Singapore, Malaysia, Thailand, Vietnam, Myanmar, Indonesia, Philippines, Bangladesh, Hong Kong, Taiwan, South Korea, China, Japan, India, UAE, Central and South Africa, Denmark, Germany, Italy, Lithuania, Russia, United Kingdom, and USA. Because of our worldwide presence, you can be sure that Esco is within your reach.

High Quality, Reliable, and Dependable. Esco products are of high quality, reliable, and dependable; assuring customers of research accuracy. Cross functional teams from Esco Production, R&D, Quality Assurance, and Senior Management, are regularly assembled to review and implement areas for improvement.

Esco Cares for Your Safety. Esco focuses on providing safety not just for your samples but also for you and the environment.

Esco Cares for Your Comfort. Building ergonomic designs and reducing noise levels of the units ensures comfort for our users.

Esco Cares for the Environment. One in every four of Esco's employees is involved in R&D and a number of them evaluate new components and/or designs to produce energy efficient equipment. Being GREEN is more than just modifying parts used to produce a new energy efficient technology, it is also embodied in the every aspect of the company.

Customer Service and Support. Our service does not stop once purchase has been done. Esco gives on-time customer service and offers end-user seminars, service training, preventive maintenance, and provides educational materials and informative videos.

As Esco takes the opportunity to respond to the world's needs, we aim not only to contribute in the advancement of scientific discoveries but also in making the world a safer, healthier, and better place to live in.

Products and Application

Laboratory Equipment

Sample Handling and Preparation

- Class I Biological Safety Cabinets
- Class II Biological Safety Cabinets
- Class II Type A2 Biological Safety Cabinets
- Class II Type B1 Biological Safety Cabinets
- Class II Type B2 Biological Safety Cabinets
- Class III Biological Safety Cabinets
- Horizontal Laminar Flow Cabinets
- Vertical Laminar Flow Cabinets
- Laboratory Animal Research Workstations
- Laboratory Centrifuges

Sample Cultivation

- CO₂ Incubators, Direct Heat Air-Jacketed
- CO₂ Incubators with Cooling System
- CO₂ Incubators with High Heat Sterilization
- Laboratory Shakers

Amplification and Detection

- Conventional Thermal Cyclers
- Microplate Shakers
- PCR Cabinets

Sample Storage & Sample Protection Solutions

- Ultra-low Temperature Freezers
- Lab Refrigerators and Freezers
- Sample Database Management Software
- Intelligent Remote Monitoring Application Protocol
- Remote Monitoring, Datalogging, Programming Software
- Wireless Monitoring System

Chemical Research

- Ducted Fume Hoods
- Ductless Fume Hoods
- Filtered Storage Cabinets
- Powder Weighing Balance Enclosure
- Exhaust Blowers
- Fume Hood Airflow Monitor

General Equipment

Laboratory Thermostatic Products

- Forced Convection Laboratory Oven
- Forced Convection Laboratory Incubator
- Natural Convection Laboratory Incubator
- Refrigerated Laboratory Incubator

Medical / IVF Equipment

Controlled Embryo Handling

- Esco Multi-Zone ART Workstation
- Esco Multi-Zone ART Workstation Class II
- AVT Anti-Vibration Table
- Semi-Closed Environment (SCE) IVF

Safe Embryo Culture

- MIRI® Multiroom Incubator
- MIRI® II Multiroom Incubator
- Mini MIRI® Humidified Incubator
- Mini MIRI® Dry Incubator
- CelCulture® CO₂ Incubator

Innovative Time-Lapse Imaging

MIRI® Time-Lapse Incubator

Accurate Quality Control

• MIRI® GA Gas and Temperature Validation Unit

Unique Consumables

CultureCoin[®]

Esco Pharma Products

Airflow Containment

- BioBooth™
- Ceiling Laminar Airflow (CLAF)
- Cytoculture® Cytotoxic Safety Cabinet
 Pharmacon™ Downflow Booth
- Esco Garment Storage Cabinet Esco Glassware Hoods
- Laminar Flow Horizontal/Vertical Trolley (LFH/VT)
- Laminar Flow Straddle Units
- Evidence Drying Cabinet

Isolation Containment

- Advanced Processing Platform Isolator (APPI)
- Aseptic Containment Isoaltor (ACTI)
- Blood Cell Labelling Isolator
- Streamline® Closed Restricted Access Barrier System (SLC-RABS)
- Containment Barrier Isolator (CBI)
- CBI-Unidirectional (CBI-U)
- CBI-Turbulent (CBI-T)
- CBI-Class III Biosafety Cabinet (CBI-III)
- CBI-Convertible Class III/Class I Biosafety Cabinet (CBI-H)
- Isoclean® Healthcare Platform Isolator (HPI)
- HPI-G3-Without Filter Below Work Zone - HPI-G3-With Filter Below Work Zone - HPI-Inflatable Seal (HPI-IS)
- · General Processing Platform Isolator
- GPPI-Inflatable Seal (GPPI-IS)
- GPPI-Static Seal (GPPI-SS)
 Streamline® Compounding Isolator
- SCI Isolator Configuration
- SCI Class III Biosafety Cabinet (SCI-III) Technetium Dispensing Isolator
- Turbulent Flow Aseptic Isolator (TFAI)
- Weighing and Dispensing Containment Isolator (WDCI)

Cross Contamination Facility Integrated Barrier

- BioPass™ Pass Through
 Cleanroom Air Showers
- Dynamic Pass Boxes/ Dynamic Floor Laminar Hatches
- Infinity® Air Shower Pass Box

- Esco Sputum Booth
 Infinity® Pass Boxes
 Infinity® Cleanroom Transfer Hatch
- Soft Capsule® Soft Wall Cleanroom

Ventilation Containment

Ventilated Balance Enclosure

Esco VacciXcell Products

Bioreactors and Fermenters

- CelXrocker[™]
- CelCradle[™]
- CelShaker[™] CelCradle™ X
- CelCradle Semi-Automated Harvesting System[™] (CCX-SAH)
- BioXcell™
- StirCradle[™]
- StirCradle[™] PRO
- TideXcell™
- TideXcell™ Harvesting System (TXLHS)
- VXL[™] Hybrid Bioreactor

Cell Culture Monitoring, Media and

Consumables

- Super Plus[™]
- Plus[™] Vero • Plus™ MDCK
- Plus[™] MDCK II
- BioNOC™ II macrocarriers
- GlucCell™ Glucose Monitoring System
- CVD Kit

Filling Line Equipment

- · Filling Line Isolators
- cRabs (close restricted access barriers)
- oRabs (open restricted access barriers)

Integrated Solutions

- Cell Processing Isolator
- Cell Processing Center

Esco TaPestle Rx Products

Pharmacy Compounding Solutions

- Compounding Pharmacy Isolators (SCI, HPI, CBI, GPPI)
- Safety Cabinets and Enclosures (CYT, Class II BSC, VBE, LFC)
- Aseptic Filling Systems

Radiopharmacy Equipment Solutions

- Radioisotope Fume Hood
- Lead-lined Biosafety Cabinet
- Technetium Dispensing Isolator
- Blood Cell Labeling Isolator • GMP-compliant Radioisotope Dispensing Isolator



CelCulture®

CO₂ Incubators with High Heat Sterilization

INTRODUCTION

Introducing Esco's CelCulture® CO_2 Incubator with 180 °C High Heat Sterilization Cycle, offering efficient contamination protection and hassle-free maintenance without compromising accuracy and reliability in maintaining optimal conditions for cell growth.

The CelCulture® CO₂ Incubator has more design configurations suitable to meet the demands of every cell culture laboratory, taking your scientific dreams a step closer to reality.

NEW FEATURES

180°C HIGH HEAT STERILIZATION

Quick and hassle-free elimination of contaminants in the chamber and its interior components.

HEAT-RESISTANT SENSORS

Maintenance-free sensors are to be included during sterilization.

TEMPERATURE FAIL-SAFE SYSTEM

Over-temperature protection device prevents overshooting of temperature to + 0.4 °C of the set point.

WATCHDOG SYSTEM-FAILURE MODE

The auto-reset watchdog will automatically reset the system in the unlikely event of system failure, preventing the controller from freezing.

%CO₂ FAILURE MODE PROTECTION

Prevents build-up of $\%CO_2$ over set point in cases of CO_2 sensor defect. The system will automatically stop the valve from injecting CO_2 after a certain period.

Available in 170 L (6.0 ft³) and 240 L (8.5 ft³) compact footprints

ULPA FILTER

- 99.999% efficient, superior to conventional HEPA filters
- Filters air continuously
- · Chamber returns to ISO Class 5 cleanliness in 11 minutes upon door closing to prevent contamination



SHELVING -

- Perforated shelving to improve uniformity
- Anti-tip
- Stainless steel
- · Built-in grip
- Dismantles without tools for easy cleaning

DIRECT HEAT & AIR JACKET

- Fast and uniform heating
- Rapid temperature recovery
- Air jacket improves chamber stability



DUCT WORK -

- Directs air flow for rapid recovery and excellent uniformity
- · Easily removed for cleaning




WATER PAN

- Precisely heated by base heater to provide high humidity
- Gentle airflow over water surface accelerates humidity recovery



ROUNDED CORNERS

- Seamless design
- Facilitates easier cleaning




• Ensures a healthier, safer, and cleaner lab environment.




IMPROVED CULTURING ATMOSPHERE FOR BETTER CELL GROWTH


VIVOCELLTM PRECISE PARAMETER CONTROL

Direct heat and air jacketed design allows even distribution of heat with less than ±0.35°C* temperature variation at 27 points in the chamber, following **DIN 12880: 2007** testing standards.






VENTIFLOW™ FORCED CONVECTION




Gentle airflow accelerates homogenization and filtration of chamber atmosphere, preventing dehydration of samples while minimizing sample stress. Blower fan automatically stops when main door is opened to minimize contamination risk.

FAST PARAMETER RECOVERY



Precise and stable sensor system combined with the SmartSense™ microcontroller allows quick parameter recovery without overshooting.

COMPLETE CONTAMINATION CONTROL




Complete Cycle lasts up to 12 hours.

8 EFFICIENCY FILTER HEP/ 0.12 MOST PENETRATING PARTICLE SIZE (MICROMETERS)

ULPA FILTRATION SYSTEM

Has 10x more filtering efficiency than HEPA filter for a cleaner and safer chamber atmosphere.




ISOCIDE™ ANTIMICROBIAL **SURFACE COATING**

Enhances sample protection by inhibiting microbial growth on the external surfaces.

180°C HIGH HEAT STERILIZATION

Conforms to the International Standards for dry heat sterilization and proven to be effective in killing normally-resistant fungi, bacterial spore, and vegetative cells. Nontoxic and noncorrosive sterilization that completes within 12 hours leaving the chamber cool and dry at the end of the cycle.

REAR PANEL





Power Supply Inlet

Connects the incubator unit to the power source.



6 N₂ Gas Supply Inlet (for Suppressed O, model)

Only applicable for models with N₂* control function. Inlet pressure requirement is 15 psi.

* O₂ and N₂ functions are applicable only to models with Suppressed O₂.



Cooling Fan

Prevents the electrical panel from overheating.



CO, Gas Supply Inlet

Connects the CO₂ gas supply to the incubator. Inlet pressure requirement is 15 psi.



RS485 Communication Port

Provides serial communication port for PC. It can be daisy-chained from one product to another and can also be connected to a PC



8 Access Ports

Allows cables, hoses or additional sensors to be routed into the work space. A rubber stopper is installed as standard configuration and is part of standard accessories.



4 Analog Port (Optional)

Allows the incubator to output analog signals representing temperature, CO₂/O₂* concentration and relative humidity, depending on the options available in the incubator. This allows the incubator to be connected to an inhouse data acquisition or alarm system.



9 0.2µm Gas Inlet Filter

Provided to remove any contaminants from the gas supply.



5 Alarm Contact

A set of relay contacts located on the rear of the unit is provided to monitor temperature, humidity or CO₂ alarms. The alarm contacts can be connected to a remote alarm system.







ENGINEERING DRAWING


Front view 660 mm (26.0") ė:∃ **;;;** ė 906 mm (35.6") MODEL 170 L 505 mm (19.9") 750 mm (29.5")


Side view





Rear view





- Control panel
 On / off switch
- Blower fan
 ULPA filter

(35.6")

906 mm (

MODEL

240 L

- 5. Sensors
- Access port
- 7. Adjustable shelves 8. Humidity pan
- 9. CO₂ gas supply 10. Alarm contact
- 11. Analog output 12. RS485

13. Cooling fan 14. Power Supply Inlet

ORDERING INFORMATION

595 mm (23.4")

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER			
MODEL ITEM CODE DESCRIPTION		DESCRIPTION	
CCL-170B-8-HHS	2170295	CelCulture® Incubator 170 L IR Sensor, CO ₂ Control, ULPA, 180°C HHS, 230 VAC 50/60 Hz	
CCL-240B-8-HHS	2170270	CelCulture® Incubator 240 L IR Sensor, CO ₂ Control, ULPA, 180°C HHS, 230 VAC 50/60 Hz	

SUPPRESSED O ₂ MODEL WITH STAINLESS STEEL CHAMBER			
MODEL	ITEM CODE	DESCRIPTION	
CCL-170T-8-HHS	2170297	CelCulture® Incubator 170L IR Sensor, CO ₂ /O ₂ Control, ULPA, 180°C HHS, 230 VAC 50/60 Hz	
CCL-240T-8-HHS	2170300	CelCulture® Incubator 240L IR Sensor, CO ₂ /O ₂ Control, ULPA, 180°C HHS, 230 VAC 50/60 Hz	

IR SENSOR MODEL WITH STAINLESS STEEL CHAMBER (NO ULPA FILTER)			
MODEL	ITEM CODE	DESCRIPTION	
CCL-170B-8-NF-HHS	2170298	CelCulture® Incubator 170 L IR Sensor, CO ₂ Control, 180°C HHS, 230 VAC 50/60 Hz, No ULPA Filter	
CCL-240B-8-NF-HHS	2170299	CelCulture® Incubator 240 L IR Sensor, CO ₂ Control, 180°C HHS, 230 VAC 50/60 Hz, No ULPA Filter	

SUPPRESSED O ₂ MODEL WITH STAINLESS STEEL CHAMBER (NO ULPA FILTER)			
MODEL	ITEM CODE	DESCRIPTION	
CCL-170T-8-NF-HHS	2170301	CelCulture® Incubator 170 L IR Sensor, CO ₂ /O ₂ Control, 180°C HHS, 230 VAC 50/60 Hz, No ULPA Filter	
CCL-240T-8-NF-HHS	2170302	CelCulture® Incubator 240 L IR Sensor, CO ₂ /O ₂ Control, 180°C HHS, 230 VAC 50/60 Hz, No ULPA Filter	



GENERAL SPECIFICATIONS		CCL-170HHS	CCL-240HHS
		TEMPERATURE	
Temperature Control Method		Direct Heat and Air Jacket using Microcontroller PI	
Ambient Temperature Range		18 to 30 °C (6	4.4 to 86.0 °F)
Temperature Co	ontrol Range, °C	Ambient	+7 to 60
Temperature U	niformity, °C *	Standard Unit: <± 0.35 Suppressed O₂ model: <± 0.4	
Temperature Fl	uctuation, °C *	±C	.2
	ecovery Time** ds door opening, 98% from initial value)	≤5 mi	nutes
		CO ₂	
CO₂ Control Sys	tem	Microcor	ntroller PI
CO₂ Control Ra	nge, % CO ₂	0.1-	19.5
CO₂ Fluctuation	ı, % CO ₂ ***	± 0.2	± 0.3
CO ₂ Sensor		Infrared (I	R) Sensor
CO₂ Recovery T (after 30 secon	ime**** ds door opening, 98% from initial value)	At 5.0% CO₂ by volume (Standard unit): ≤5 minutes Suppressed O₂ model: ≤8 minutes	At 5.0% CO₂ by volume (Standard unit): ≤5 minutes Suppressed O₂ model: ≤8 minutes
		O ₂	
O ₂ Control Syst		Microcor	ntroller PI
O ₂ Control Range, % O ₂		1-	18
O ₂ Fluctuation,	% O ₂ ****	± ().3
O₂ Sensor		Zirconia O₂ Sensor	
O, Recovery Time****** (after 30 seconds door opening, 98% from initial value)		At 5.0% O₂ by volume: ≤12 minutes	At 5.0% O₂ by volume: ≤12 minutes
		HUMIDITY	
Humidification	Method	Humidity pan	
Humidity Rang	e (at 37°C)******	85-90	
		PHYSICAL CONSTRUCTION	
Interior Volume		170 L (6 ft³) 240 L (8.5 ft³)	
External Dimen	sions (W x D x H)	660 x 660 x 906 mm (26.0" x 26.0" x 35.6")	750 x 770 x 906 mm (29.5" x 30.3" x 35.6")
Internal Dimen	sions (W x D x H)	505 x 535 x 633 mm (19.9" x 21.1" x 24.9")	595 x 640 x 633 mm (23.4" x 25.2" x 24.9")
Net Weight		101 kg (222.7 lbs.)	121 kg (266.8 lbs.)
	Main Body	Electrogalvanized steel with Is	ocide™ antimicrobial coating
	Interior Material	Stainless ste	el, type 304
Chamber	Number of Shelves	4	1
Construction	Maximum Number of Shelves	7	7
	Shelves Area (W x D)	465 x 470 mm (18.3" x 18.5")	550 x 560 mm (21.7" x 22.0")
	Maximum Load per Shelf	11 kg/shelf (24.3 lbs./shelf)	15 kg/shelf (33.1 lbs./shelf)
Electrical	Nominal Power at 37°C	42.2 W	42.2 W
Configuration 220-240 VAC,	Maximum Power Consumption	1300 W	1500 W
50/60 Hz	Full Load Amps	5 A	7 A
Shipping Weight		140 kg (308.6 lbs)	160 kg (352.7 lbs)
Shipping Dimensions (W x D x H)		850 x 720 x 1120 mm (33.5" x 28.3" x 44.1")	850 x 850 x 1120 mm (33.5" x 33.5" x 44.1")
Shipping Volun	ne	0.70 m³ (24.85 ft³)	0.79 m³ (28.03 ft³)
		CONTAMINATION CONTROL	
Contamination Control Methods		 Main body is electrogalvanized steel with Isocide™ antimicrobial coating; 180°C high heat sterilization cycle; ULPA filter (optional) - filter must be removed during decon 0.2 µm gas inlet filter 1-micron air circulation filter 	

OPTIONS AND ACCESSORIES

	DESCRIPTION	COA CODE	ITEM CODE
	HUMIDITY DISPLAY This option allows the incubator to monitor the relative humidity inside the chamber. The sensor is easy to install and has excellent accuracy. The	COA-1001 (factory-installed)	5170470
	airflow in the chamber does not affect the measurement. The sensor is maintenance-free and does not need to be removed prior to sterilization.	COA-1001-F (field-installed)	5170471
0	CO ₂ BACKUP This option allows two tanks of CO ₂ to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low	COA-1002 (factory-installed)	5170472
	gas pressure is detected on the primary tank.	COA-1002-F (field-installed)	5170473
	N_2 BACKUP This option allows two tanks of N_2 to be connected to the incubator. It will automatically switch from the primary tank to the secondary tank when low	COA-1007 (factory-installed)	5170490
	gas pressure is detected on the primary tank.	COA-1007-F (field-installed)	5170491
	ANALOG OUTPUT A set of relay contacts is provided at the rear of the incubator that allows the incubator to output analog signals representing the temperature, $\%CO_2$, $\%O_2$ and relative humidity, depending on the options available in the incubator. This allows the chamber to be connected to an in-house data acquisition or alarm system. This option can also be field-installed.	COA-1005 (factory-installed)	5170475
	The analog signal outputs can be set to operate in either voltage DC (0-5 VDC) or current (4-20 mA) mode. The factory default setting is voltage. Switch on the board to toggle between the modes.	COA-1005-F (field-installed)	5170476
	2-STAGE GAS REGULATOR FOR CO_2/N_2 CO_2 and N_2 gas input regulators reduce pressure from the tank to the incubator. It has dual pressure gauges, barbed line connection and shutoff valve. It prevents over-pressurization of the gas supply into the incubator which could cause the tubing to burst.	COA-2005-F	5170481
	EXTRA STAINLESS STEEL SHELF Each CO_2 incubator comes standard with 4 shelves and it can accommodate up to a maximum of 7 shelves.	COA-2007-F (for 170 L models)	5170327
		COA-2025-F (for 240 L models)	5170426
	ROLLER BASE Roller base is available with casters for mobility of your incubators and to provide protection against floor contamination.	COA-2001-F (for 170 L models)	5170478
		COA-2019-F (for 240 L models)	5170420
	FLOOR STAND 200 MM (8.0") WITH ADJUSTABLE FEET Floor stands are available with adjustable feet, with a nominal range of 180 mm to 250 mm (7.1" to 9.8") for comfortable access to the incubator and to avoid floor contamination.	COA-2002-F (for 170 L models)	5170479
	to avoid 11001 Contamination.	COA-2021-F (for 240 L models)	5170422
	FLOOR STAND 700 MM (27.6") WITH CASTERS This support stand raises the incubator to a height of 700 mm (27.6") above the floor for comfortable access. It comes with casters for mobility of your incubators.	COA-2003-F (for 170 L models)	5170480
		COA-2023-F (for 240 L models)	5170424

	DESCRIPTION	COA CODE	ITEM CODE
	2-UNITS FLOOR STAND STACKING KIT This floor stand allows two incubator units to be stacked without being physically in contact with each other. For the lower unit, it uses roller base for mobility and for easy pull out of the lower unit in case of troubleshooting. Floor stand for upper unit also has casters for easy relocation.	COA-2004-F (for 170 L models)	5170489
		COA-2042-F (for 240 L models)	5170999
	ELECTRONIC CO ₂ ANALYZER, FOR CO ₂ / TEMP MEASUREMENT (WITH TEMPERATURE PROBE)	COA-2010-F	5170329
	ELECTRONIC CO ₂ + O ₂ ANALYZER, FOR CO ₂ / O ₂ / TEMP MEASUREMENT (WITH TEMPERATURE PROBE)	COA-2016-F	5170397
	ELECTRONIC CO ₂ + O ₂ + RH ANALYZER, FOR CO ₂ / O ₂ / RH / TEMP MEASUREMENT (WITH TEMPERATURE PROBE)	COA-2017-F	5170398
	6" CHART RECORDER, TEMP, 115/230 VAC, 50/60 HZ The chart recorder provides an easy-to-read graph of data vs. time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature data.	COA-2012-F	2170021
	8" CHART RECORDER, TEMP/TEMP, 115/230 VAC, 50/60 HZ The chart recorder provides an easy-to-read graph of data vs. time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 8" chart of temperature data and comes with 2 remote probes for dual temperature monitoring.	COA-2013-F	2170022
	6" CHART RECORDER, TEMP/RH, 115/230 VAC, 50/60 HZ The chart recorder provides an easy-to-read graph of data vs. time. It is a reliable, accurate, and stable instrument for on-the-spot written documentation of incubator chamber temperature. This model offers 6" chart of temperature and humidity data.	COA-2014-F	2170023
	REVERSED DOOR SWING The incubator has a door opening on the left side by default. This option allows the doors to be factory-installed as opening from the right side.	COA-1004 (factory-installed)	5170474
FOR THE STATE OF T	IQ / OQ DOCUMENTATION The execution of the IQ / OQ verifies that the incubator is installed and is operating pursuant to the validated Standard Operating Procedures (SOPs).	COA-2011-F	2170020
Voyages	VOYAGER SOFTWARE KIT Esco Voyager is a PC-based software package developed for the remote monitoring, data logging, and programming / device configuration of Esco controlled environment laboratory equipment. Compatible equipment includes laboratory ovens and incubators, low temperature incubators, CO ₂ incubators, and ultra-low temperature freezers.	Voyager	5250001

ESCO LIFESCIENCES GROUP

42 LOCATIONS IN 21 COUNTRIES ALL OVER THE WORLD



Follow us on social media, download our apps, and scan the OR code for more info.



























Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escolifesciences.com www.escolifesciences.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA Tel: +1 215-441-9661 • Fax 484-698-7757 eti.admin@escolifesciences.com

Esco Lifesciences Group Offices: Bangladesh | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Myanmar | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | UK | USA | Vietnam

